

Electron Spin Resonance

Nathan Grouse

Lisa Tran


May 17, 2011

1 Introduction

An oscillating magnetic field is used to induce transitions between energy levels taken on by electrons in this field.

1.1 Apparatus

Helmholtz coils wired in parallel, a multimeter to measure the current in the coils, a power source provides a DC voltage with a possible modulation voltage, and an oscilloscope are used.



2 Theory

Spin angular momentum:

$$\vec{S} = \sqrt{s(s+1)} \frac{h}{2\pi},$$

where $s = \frac{1}{2}$ is the spin quantum number. This makes sense considering the half integer multiples of Planck's constant.

The spin angular momentum is related to the magnetic moment:

$$\vec{\mu} = -(g_s \mu_B \vec{S}) / \frac{h}{2\pi}$$

Meaning magnetic moments also come in the same multiples as the spin angular momentum. g_s is the g-factor and its value depends on the particle being described. $g_{proton} = 5.586$, $g_{electron} = 2.0023$, and $g_{neutron} = -3.826$.

In a magnetic field that points in one direction, apparently the spin angular momentum can only have projections onto the respective axis, resulting in two values for the spin magnetic moment:

$$S_z = \hbar m_s,$$

where $m_s = \pm \frac{1}{2}$ are the possible values for the magnetic spin quantum number. So the spin magnetic moment also takes on two different values in one dimensional magnetic fields:

$$\mu_z = -(g_s \mu_B S_z) / \hbar = \pm \frac{1}{2} g_s \mu_B.$$


Two possible energy levels are defined by this too:

$$U = -\vec{\mu} \cdot \vec{B} = \pm \frac{1}{2} g_s \mu_B B,$$

$$\Delta E = h\nu = g_s \mu_B B.$$

Where ν is the minimum frequency of radiation that would induce or be emitted during a change of state for a particle in one dimension, derived from the Einstein relation.

3 Data



4 Calculations

Measuring g :

$$\nu =$$

$$g \mu_B \bar{h} B$$

$$\frac{g\mu_B}{h}(10^9) = 13.66 \frac{MHz}{mT} = 1.366 \times 10^{10} \frac{Hz}{T}$$

$$g = (1.366 \times 10^{10}) \frac{h}{\mu_B} = .975$$

Measuring δB :

$$\delta W = 1.2 cm$$

$$I_{mod,rms} = .33 A$$

$$I_{mod,p-p} = 2\sqrt{2}I_{mod,rms} = .93 A$$

$$\delta(2I_0) = I_{mod,p-p} \frac{\delta W}{10} = .011 A$$

$$\delta B = 2.115 \delta(2I_0) = .0237 mT$$

5 Error Analysis

The value of δB is within the designated range of line widths quoted in the literature. I won't check it's percent error but .05 MHz, as indicated by the decimal accuracy of the machine's readings. The uncertainty in the current readings is $\pm .005$ A, as indicated by the decimal accuracy of the multimeter. Uncertainty propagated linearly through calculating values of the magnetic field between the coils from current:

$$B_{max} = 2.115(2(I_0 + .005A)) = B + .02115mT$$

$$B_{min} = 2.115(2(I_0 - .005A)) = B - .02115mT$$

$$\delta B = \pm .02115mT$$

Percent Error on g:

$$\frac{|2.0023 - .975|}{2.0023} = 51.3\%$$

6 Conclusion

The numbers I obtained for both g and δB seem reasonable and consistent with the accepted value. I did see what I expected to see in terms of the resonance curves.

7 Questions

1. The manufacturer designed this experiment with the coils connected in parallel. A series connection would be better. Why?

If the inductances of two coils located in each others magnetic field are not equal, this can cause short circuiting.

2. The p-p modulation current $\delta(2I_0)$ for the half-width δB is obtained from:

$$\delta(2I_0) = I_{mod,p-p} \frac{\delta W}{10}$$

Where does the divisor 10 come from?

It changes the value of δW so that it represents meters.

3. In the method given for measuring δB , the scope controls are not used in a calibrated mode. Why is this OK?

We account for it by dividing δW by 10.

4. Why is the multimeter set for DC amperes for measuring g and for AC amperes for measuring the line width?

Not sure. Maybe because AC amperes are required for generating p-p values.

5. Is there an RF electric field associated with the RF coil? If so, make a sketch of what the fields look like.

I think there is an RF electric field but I don't know what it looks like.